




















Connected Care: Continuous Glucose Monitor Guide














<u>Continuous Glucose Monitor</u>	<u>Compatible Smart Device</u>	<u>Mobile App</u>
<p><u>Dexcom G6</u></p>  <p>10-day sensor 90-day transmitter 2-hour warm-up time Indicated for age 2+</p>	 <p>Visit https://Dexcom.com/compatibility for a full list of compatible smart devices</p>	 <p>Dexcom G6 Dexcom Follow</p>  <p>Dexcom Clarity Dexcom Follow (Enhanced)</p>
<p><u>Dexcom G7 10-Day Sensor</u></p>  <p>10-day sensor w/built-in transmitter 30-minute warm-up time 12-hour grace period Indicated for age 2+</p>	 <p>Visit https://Dexcom.com/compatibility for a full list of compatible smart devices</p>	 <p>Dexcom G7 Dexcom Follow</p>  <p>Dexcom Clarity Dexcom Follow (Enhanced)</p>







Connected Care: Continuous Glucose Monitor Guide

Continuous Glucose Monitor	Compatible Smart Device	Mobile App
<p><u>Dexcom G7 15-Day Sensor</u></p>  <p>15-day sensor w/built-in transmitter 30-minute warm-up time 12-hour grace period *Indicated for 18+ Coming 2nd half of 2025</p>	 <p>Visit https://Dexcom.com/compatibility for a full list of compatible smart devices</p>	 <p>Dexcom G7</p>  <p>Dexcom Follow</p>  <p>Dexcom Clarity</p>  <p>Dexcom Follow (Enhanced)</p>
<p><u>Freestyle Libre 2 Plus</u></p>  <p>15-day sensor w/built-in transmitter 60-minute warm-up time Indicated for age 2+</p>	 <p>Visit https://FreestyleLibre.us for a full list of compatible smart devices</p>	 <p>Freestyle Libre 2 (Sensor is flash scanner)</p>  <p>LibreLinkUp</p>  <p>Libre by Abbott (Sensor is Real Time CGM)</p> <p>LibreView</p> <p>Website only www.libreview.com</p>







Connected Care: Continuous Glucose Monitor Guide

Continuous Glucose Monitor	Compatible Smart Device	Mobile App
<p><u>Freestyle Libre 3 Plus</u></p>  <p>15-day sensor w/built-in transmitter 60-minute warm-up time Indicated for age 2+</p>	  <p>Visit https://FreestyleLibre.us for a full list of compatible smart devices</p>	 <p>Freestyle Libre 3</p>  <p>LibreLinkUp</p>  <p>Libre by Abbott</p> <p><i>LibreView</i></p> <p>Website only www.libreview.com</p>
<p><u>Guardian Connect 4</u></p>  <p>7-day sensor w/rechargeable transmitter 2-hour warm-up time Indicated for age 7+ Integrates w/Medtronic MiniMed 780G and InPen *Not a stand-alone sensor</p>	   <p>Visit https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility for a full list of compatible smart devices</p>	 <p>Medtronic</p> <p>MiniMed Mobile US</p>

Connected Care: Continuous Glucose Monitor Guide

<u>Continuous Glucose Monitor</u>	<u>Compatible Smart Device</u>	<u>Mobile App</u>
<p><u>Simplera</u></p>  <p>6-day sensor w/24-hour grace period 2-hour warm-up time Integrates with InPen Indicated for ages 7+ Can be worn as stand-alone sensor</p>	 <p>Visit https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility for a full list of compatible smart devices</p>	 <p>Simplera CGM US InPen Smart Pen</p>
<p><u>Simplera Sync</u></p>  <p>6-day sensor w/24-hour grace per 2-hour warm-up time Integrates with MiniMed 780G Indicated for ages 7+ *Not a stand-alone sensor</p>	 <p>Visit https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility for a full list of compatible smart devices</p>	 <p>MiniMed Mobile US</p>

Connected Care: Continuous Glucose Monitor Guide

Continuous Glucose Monitor	Compatible Smart Device	Mobile App
<p><u>Instinct by Abbott</u></p>  <p>15-day sensor w/built in transmitter 1-hour warm-up time Compatible with MiniMed 780G Indicated for ages 7+ *Not a stand-alone sensor</p>	 <p>Visit https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility for a full list of compatible smart devices</p>	 <p>MiniMed Mobile US</p>
<p><u>Senseonics Eversense 365</u></p>  <p>Sensor implanted under the skin and lasts 1 year Transmitter worn over and powers sensor</p> <ul style="list-style-type: none"> Rechargeable: 15 minutes = 24 hours Removable <p>Requires calibration with fingerstick 24-hour warm-up phase Will integrate with Sequel twist insulin pump</p>	 <p>Visit https://global.eversenseddiabetes.com/compatibility for a full list of compatible smart devices</p>	 <p>Eversense 365</p>