

Connected Care: Infusion Sets

<u>Tandem</u>	
<u>AutoSoft 90</u> 	• 6- or 9- mm flexible soft cannula • 90-degree insertion angle • 23- or 43-inch tubing • Automatic insertion device • Quick reversible connector • Built-in needle cover • Available in three colors: blue, pink, and gray • Change frequency: 48 hours if using Humalog, 72 hours if using Novolog
<u>AutoSoft 30</u> 	• 13 mm flexible soft cannula • 30-degree insertion angle • 23- or 43-inch tubing length • Automatic insertion device • Quick reversible connector • Large grip for easy disconnection • Large adhesive for extra security • Window to see insertion point • Change frequency: 48 hours if using Humalog, 72 hours if using Novolog
<u>AutoSoft XC</u> 	• 6- or 9- mm flexible soft cannula • 90-degree insertion angle • 5, 23, 32, or 43-inch tubing length • Automatic insertion device • Quick reversible connector • Easy disconnection • Built-in needle cover • Available in three colors: blue, pink, and gray • Change frequency: 48 hours if using Humalog, 72 hours if using Novolog *5-inch tubing option is only available in 6 mm cannula length, one color (gray) and only compatible with Tandem Mobi AID System*
<u>VariSoft</u> 	• 13- or 17- mm flexible soft cannula • Variable insertion angle: 20°-45° • 23-, 32-, or 43-inch tubing length • Manual insertion • Audible connector • Large adhesive for extra security • Slim profile for discreetness • Window to see insertion point • Change frequency: 48 hours if using Humalog, 72 hours if using Novolog

Connected Care: Infusion Sets

Tandem	
<u>TruSteel</u> 	• 6- or 8-mm needle • 29-gauge stainless steel needle • 90-degree insertion angle • 23- or 32-inch tubing length • Manual insertion • Large grip for easy disconnection • Second pad for extra security • Change frequency: 24 - 48 hours
<u>SteadySet</u> 	• Pre-loaded/All-In-One Inserter • 3-step insertion process • Up to 7 days of wear • Soft & kink-resistant cannula • 35° insertion angle • Ability to temporarily disconnect • Multiple tubing lengths • Available early 2026
<u>SteadySet</u> 	• Tubeless (for Tandem Mobi) • Independent cartridge/infusion set • Ability to temporarily disconnect • Pre-loaded/All-In-One Inserter • 3-step insertion process • Up to 7 days of wear • Soft & kink-resistant cannula • 35° insertion angle • Available early 2026


Connected Care: Infusion Sets

<u>Medtronic</u>	
<u>Medtronic Extended infusion set</u> 	• All-in-one insertion device • 25 guage soft cannula • 90° insertion angle • 6- or 9- mm cannula length • 23- or 32-inch tubing length • Change frequency: up to 7 days
<u>MiniMed Mio Advance</u> 	• All-in-one insertion device • 25 guage soft cannula • 90° insertion angle • 6- or 9- mm cannula length • 23- or 43-inch tubing length • Change frequency: 2 - 3 days
<u>MiniMed Quick-set</u> 	• Quick-serter or Manual insertion • 25 guage soft cannula • 90° insertion angle • 6- or 9- mm cannula length • 18, 23, 32, or 43-inch tubing length • Change frequency: 2 - 3 days
<u>MiniMed Silhouette</u> 	• MiniMed Sil-serter or Manual insertion • 25 guage soft cannula • 30° - 45° insertion angle • 13- or 17- mm cannula length • 18, 23, 32, or 43-inch tubing length • Change frequency: 2 - 3 days
<u>MiniMed Sure-T</u> 	• Manual insertion • 29 guage steel needle/cannula • 6-, 8-, or 10-mm needle/cannula • 90-degree insertion angle • 18, 23, or 32-inch tubing length • Second pad for extra security • Change frequency: 24 - 48 hours

Connected Care: Infusion Sets

Omnipod

Omnipod 5



- Pod contains cannula and reservoir
 - Needle places cannula under skin then retracts into pod
- 45° insertion angle
- 6 mm soft cannula
- No tubing
- Waterproof up to 25 feet for 60 minutes
- Change frequency: Up to 72 hours
 - 8-hour grace period with Omnipod 5 if insulin available in pod

Connected Care: Infusion Sets

iLet

Inset I



- Automatic Insertion Device
- 6 mm, 25-gauge Teflon cannula
- 90° insertion angle
- 23- & 32-inch tubing length
- Change frequency: 2 – 3 days

Contact Detach



- Manual Insertion
- 6 mm, 29-gauge, Stainless Steel cannula
- 90° insertion angle
- 23-inch tubing length
- Change frequency: 1 - 2 days

Connected Care: Infusion Sets

twiist

Cleo 90



- All-in-one insertion device
- 360° on-site rotation for increased flexibility
- 25 guage soft cannula
- 90° insertion angle
- 6 mm cannula length
- 24-inch tubing length
- Change frequency: 2 -3 days