

## Connected Care

### Quick Guide to Starting a Continuous Glucose Monitor

Please download the following apps if using smartphone as receiver



#### Dexcom G6



- Allows users to see their glucose readings on smartphone/device
- Allow users to share data with up to ten followers
- Able to use app, receiver, and AID system at the same time
- Download to user's smartphone/device

#### Dexcom G7



- Allows users to see their glucose readings on smartphones/device
- Allow users to share data with up to ten followers
- Able to use app, receiver, and AID system at the same time
- Download to user's smartphone/device

#### Dexcom Clarity



- User/parent/caregiver (over 18) must create account in their name
  - If user is a minor child, parent/caregiver must create dependant account
- Allows user to generate reports to view blood glucose data and share with provider
- Download to user/parent/caregiver smartphone/device
- **Clinic Clarity Share Code:** \_\_\_\_\_

#### Dexcom Follow



- Allows parent/caregiver to follow user's CGM data
- Available when using a smartphone/device as the receiver
- Download to parent/caregiver smartphone/device

#### Dexcom Follow (Enhanced)



- Allows parent/guardian to follow user's CGM data
- Available when using smartphone/device as the receiver
- Allows user to see data from twenty-five Sharers
  - User can notate relationship to Sharer

Additional training videos can be found at:

<https://www.dexcom.com/training-videos>

For a list of compatible smart devices please visit [Dexcom.com/compatibility](https://www.dexcom.com/compatibility)

## Connected Care

### Quick Guide to Starting a Continuous Glucose Monitor

Please download the following apps if using smartphone as receiver



#### Freestyle Libre 2



- Allows users to see their glucose readings on smartphone/device
  - If using Libre 2 app, CGM must be scanned once every 8 hours
- Download to user's smartphone/device
- Cannot use app, receiver, and AID system at the same time
- Allows users to share data with up to twenty followers

#### Freestyle Libre 3



- Allows users to see their glucose readings on smartphone/device
- Allows users to share data with up to twenty followers
- Cannot use app, receiver, and AID system at the same time
- Download to user's smartphone/device

#### Libre by Abbott



- Allows user to see glucose readings on smart phone/device
- Compatible with Freestyle Libre 2 Plus and Freestyle Libre 3 Plus
  - Turns Freestyle Libre 2 Plus into real time CGM (no scanning required)
  - Allows for more alarm customization than Libre 2 and Libre 3 apps
- Cannot use app, receiver, and AID system at the same time
- Download to user's smart phone/device

#### LibreLinkUp



- Allows parent/guardian to follow user's CGM data
  - Must create LibreLinkUp account
- Download to user or parent/guardian smart phone/device

#### LibreView

LibreView

- Allows user to generate reports to view blood glucose data and to share with provider [www.libreview.com](http://www.libreview.com) (only available as website)
- **Clinic Share Code:** \_\_\_\_\_

Additional training videos can be found at:

<https://freestyleprovider.abbott/us-en/freestyle-libre-resources.html>

For a list of compatible smart devices please visit [FreestyleLibre.us](http://FreestyleLibre.us)