












Connected Care: Continuous Glucose Monitor Guide











Continuous Glucose Monitor	Insulin Delivery Device	Smart Device	Mobile Apps
<p><u>Dexcom G6</u></p>  <p>10-day sensor 90-day transmitter 2-hour warm-up time Indicated for age 2+</p>		 	
<p><u>Dexcom G7 10-Day Sensor</u></p>  <p>10-day sensor w/built-in transmitter 30-minute warm-up time 12-hour grace period Indicated for age 2+</p>		  	

Visit <https://Dexcom.com/compatibility> for a full list of compatible smart devices

Connected Care: Continuous Glucose Monitor Guide

















Continuous Glucose Monitor	Insulin Delivery Device	Smart Device	Mobile Apps
<p><u>Dexcom G7 15-Day Sensor</u></p>  <p>15-day sensor w/built-in transmitter 12-hour grace period 60-minute warm-up time Indicated for 18+</p>			
<p><u>Freestyle Libre 2 Plus</u></p>  <p>15-day sensor w/built-in transmitter 60-minute warm-up time Indicated for age 2+</p>			 <p>Flash-Scanner CGM Real-Time CGM</p> <p>Website only www.libreview.com</p>

Visit <https://Dexcom.com/compatibility> for a full list of compatible smart devices

Visit <https://FreestyleLibre.us> for a full list of compatible smart devices

Connected Care: Continuous Glucose Monitor Guide











Continuous Glucose Monitor	Insulin Delivery Device	Smart Device	Mobile Apps
<p><u>Freestyle Libre 3 Plus</u></p>  <p>15-day sensor w/built-in transmitter 60-minute warm-up time Indicated for age 2+</p>		 	   <p>LibreView</p> <p>Website only www.libreview.com</p>
<p><u>Guardian Connect 4</u></p>  <p>7-day sensor w/rechargeable transmitter 2-hour warm-up time Indicated for age 7+ *Not a stand-alone sensor</p>		  	 

Visit <https://FreestyleLibre.us> for a full list of compatible smart devices

Visit <https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility> for a full list of compatible smart devices

Connected Care: Continuous Glucose Monitor Guide











Continuous Glucose Monitor	Insulin Delivery Device	Smart Device	Mobile Apps
<p><u>Simplera</u></p>  <p>6-day sensor w/24-hour grace period 2-hour warm-up time Indicated for ages 7+ Can be worn as stand-alone sensor</p>			
<p><u>Simplera Sync</u></p>  <p>6-day sensor w/24-hour grace period 2-hour warm-up time Indicated for ages 7+ *Not a stand-alone sensor</p>			

Visit <https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility> for a full list of compatible smart devices

Connected Care: Continuous Glucose Monitor Guide



Continuous Glucose Monitor	Insulin Delivery Device	Smart Device	Mobile Apps
<p><u>Instinct by Abbott</u></p>  <p>15-day sensor w/built in transmitter 1-hour warm-up time Indicated for ages 7+ *Not a stand-alone sensor</p>			
<p><u>Senseonics Eversense 365</u></p>  <p>Sensor implanted under the skin and lasts 1 year Transmitter worn over and powers sensor</p> <ul style="list-style-type: none"> Rechargeable – 15 minutes = 24 hours Removable <p>Requires calibration with fingerstick 24-hour warm-up phase Will integrate with Sequel twiist insulin pump</p>	 <p>Coming 2026</p>		

Visit <https://www.medtronicdiabetes.com/customer-support/app-support/device-compatibility> for a full list of compatible smart devices
 Visit <https://www.eversensecg.com/compatibility> for a full list of compatible smart devices